

CareTec IV TENS Unit

The **CareTec IV** combines the waveforms of Russian, EMS, TENS and IF therapy in a portable unit. The device offers preset modes that can be user defined. Maximum portability and simple operation offering you all the options of a clinical device in a portable package.

Features:

- Dual channel device
- Simple to operate, portable device for pain management
- Four waveforms, Russian, EMS, TENS and IF
- Comes complete with carrying case, lead wires, and 9-volt battery
- AC adapter, manual, and pack of four self adhesive reusable electrodes
- Power Source: Standard 9-volt battery or AC adapter
- Dimensions: 4.5" x 2.55" x 0.9"
- Weight: 0.28 lb. with battery
- One year warranty



	Russian	EMS	IF	TENS
Wave Form	Biphasic Square	Monophasic Square	Biphasic Square	Monophasic Square
Pulse Amplitude	0-100mA	0-100mA	0-70mA	0-60mA
Pulse Rate	2500 Hz, adjustable burst frequency	0-150 Hz	4000-4150 Hz	0-150 Hz
Pulse Width	150 μ S	50-300 μ S	125 μ S	50-300 μ S
Modes (User Defined)	3	3	4	4
Timer	1-60 & Continuous	1-60 & Continuous	1-60 & Continuous	1-60 & Continuous