

InTENSity™ 10

The InTENSity™ 10 digital TENS device allows patients to select from ten pre-set body icons for targeted therapy, making it the easiest and most effective TENS device available on the market today. With the press of two buttons, the easy-to-use device is ready to operate effectively based on the area of the body that is being treated. Comes complete with an AC adapter for effective therapy on-the-go.

Features:

- Simple to operate, portable device for pain management
- 2-button operation
- Device has ten pre-set body part selections, making it easy for patients to target certain areas in pain
- Body part selections include: neck, shoulder, mid back, low back, elbow, wrist, hip, knee, ankle, and foot
- Individual body selections have a customized treatment program based on researched protocols
- Comes complete with carrying case, two sets of lead wires, 9-volt battery, AC adapter, manual, and self-adhesive electrodes
- Power Source: Standard 9-volt battery or AC adapter
- Dimensions: 4.5" x 2.55" x 0.9"
- Weight: 0.28 lb with battery
- One year warranty



	Neck	Low Back	Mid Back	Shoulder	Elbow	Wrist	Knee	Hip	Foot	Ankle
Pulse Width	M100-150	F260	F100	F260	F100	F100	M100-150	F200	M100-160	F100
Pulse Rate	M60-100	M50-80	M100-150	M80-100	F100	F260	F120	M100-150	M60-100	F100