

MICRO 400

MENS/TENS Unit

Features:

- *Dual channel (4-lead) analog microcurrent*
- *Selectable frequencies (.3, 9.1, 30, or 292 Hz)*
- *Choice of micro-ampere or milli-ampere (TENS) currents*
- *Square or sloped waveform*
- *Adjustable amplitude*
- *Bipolar [(+) and (-)] direct current waveform*
- *LED output indicator light*
- *Durable metal belt clip*
- *1-year warranty*



Technical Specifications:

Channel:	Dual, isolated between channels
Amplitude:	Adjustable, 0 - 8 mA
Waveform:	Square or slope
Carrier Frequency:	15000 Hz
Modulating Frequency:	0.3Hz, 9.1Hz, 30Hz, 292Hz
Power Supply:	9-volt alkaline battery
Size:	2.7 x 6.2 x 9.5 cm
Weight:	150 grams including battery

Comes with 1 pkg of 4 reusable self-adhering electrodes, 1 set 43 pin lead wires, 1 each 9V heavy duty battery, carrying case, and instruction manual.

TENS  **R**
Your Real Time Healthcare Solution

Phone: 877-372-4441 | Fax: 877-372-5551